
Crizal Sapphire lenses with Scotchgard Protector – Launched April 27, 2010

Features and Benefits:

- **The Enhanced Light Transmission (ELT) System™** - An Essilor exclusive technology to further reduce glare and provide maximum light transmission through the lens. As a result, patients will benefit from **50% less reflection** compared to Crizal Avancé™ lenses with Scotchgard™ Protector.
- **Superior Crizal® Technology** – Crizal Sapphire lenses with Scotchgard Protector are powered by technology incorporating all of the best features and benefits from Crizal Avancé with Scotchgard Protector including:
 - **The SR Booster™ Layer:** The most significant breakthrough in **scratch resistance**
 - **The High Surface Density™ (HSD) Process:** The best **super-hydrophobic top coat** on the market
 - **Anti-static Properties:** A superior anti-static layer which repels dust and particles from the surface of the lens
 - **Pad Control System™:** An Essilor exclusive layer ensuring edging precision and protection

NEW Pad Control System:

This product has a slight blue reflective hue. The Pad Control System, which is normally blue in color, has been altered for this product only to an **ash/gray** color to distinguish it from the lens and must be wiped off **AFTER** edging and mounting.

Availability:

Find the availability chart within the Resource Library section of this website under “Product Availability”.