



Introducing  
EOS Wrap™  
from KBco



THE POLARIZED LENS COMPANY

Vision shouldn't be an obstacle

## Imagine-Great Optics for Oversized Fashion

And Performance Wrap Frames. KBco, The Polarized Lens Company has done it again. The polarized EOS Wrap™ PAL is a comprehensive solution. It's unique use of asphericity delivers outstanding optics and an expansive distance viewing zone in oversized wrap frames for Rx's ranging from -4.50D to +2.50D. And, only EOS Wrap comes in a blank size large enough—83mm with 3 mm of decentration—to take on oversized wrap frames. EOS Wrap is available in KBco True Grey, KBco Grey C and KBco Brown C and features our Super WearGard coating. Best of all, EOS Wrap is available from all full service laboratories.



For further information, please call (800) 722-8776

### Polycarbonate EOS Wrap Progressive

dia. (mm)	front base curve	true front curve	front curve radius (mm)	back curve	back curve radius (mm)	distance from polarizing layer to front surface (mm)	sag value 50mm bell	center thickness (mm)	edge thickness (mm)
83	8.25	8.242	64.3	8.60	61.62	0.3	5.06	13.09 +/- 0.5	13.2 +/- 0.5

\*Please do not sag, use printed sag values for power calculation as the product is an aspheric design

### Specifications

Material: Polycarbonate  
Abbe Value: 32  
Refractive Index: 1.586  
Specific Gravity: 1.20  
Minimum Center Thickness: 1.6 mm  
Color: KBco Grey-C, KBco Brown-C, KBco True Grey